

SPECIFICATIONS

Clean Flow 3003 Specifications

Total Horsepower - 3HP Triple Motor

Total air delivery at 80 -100 PSI - 15.3 CFM

Air Tank Capacity - 30 Gallons

Sound level - 62db

Maximum Users - 6 - 8

Weight and Dimensions

Weight - 298 lbs. / 340 lbs. w/pallet

Height - 33 Inches

Width - 24 Inches

Length - 41 Inches

Electrical Requirements

Voltage requirements: 240V*

3 horsepower - Triple

*240 V / 30 amp dedicated circuit

Current Amperage 7A 50/60hz

208 Line Voltage may require Buck Booster

Recommended:

Wire Size - 10 AWG

Circuit Breaker - 30 Amp



Manufactured by: Dansereau Health Products, Inc. / 1581 Commerce Street, Corona, CA 92880

Ph: (800) 423 - 5657 / Fax: (951) 549 - 1411

Website DHPDENTAL.COM